



# More reasons to trust

## Avalon FM20 and FM30 with Smart Pulse fetal and maternal monitoring

The Avalon FM series fetal and maternal monitors are Philips first to offer cross-channel verification using Smart Pulse. This innovative feature allows automated maternal pulse detection via the Toco MP transducer without the need to monitor maternal SpO<sub>2</sub> or ECG separately. Built-in NST Trace Interpretation, a clinical decision support (CDS) application for antepartum monitoring, can help to improve the efficiency of your daily work by sorting out traces that may be of concern. With the Avalon FM series, you have the confidence of continuous data with a backup memory, battery operation, and LAN interface, combined with smart transducers and a color touchscreen.

### Key advantages

- Automated coincidence detection with Smart Pulse
- Optional battery allows continuous monitoring during maternal transport within healthcare facilities
- The Avalon FM series can monitor triplets on a single monitor

**PHILIPS**  
sense and simplicity

# Smart choices at every point

Cross-channel verification continuously compares all fetal and maternal heart rates allowing each to be measured separately, enhancing confidence. Easily view maternal and fetal direct ECG waves on the large, colorful display.

These antepartum and intrapartum monitors provide:

- Separate maternal pulse measurement
- Continuous monitoring during transport in healthcare facilities
- Integrated monitoring of maternal pulse rate and blood pressure (optional)
- External monitoring of multiple fetal heart rates, uterine activity, and fetal movement
- An extensive set of internal fetal parameters, such as direct fetal heart rate and intrauterine pressure
- Optional maternal SpO<sub>2</sub>

## From the very first visit

Avalon FM fetal monitors connect easily to OB TraceVue, the Philips obstetrical surveillance and information management system.

## Enlightened enough to handle triplets on a single monitor

**LAN interface** for compatibility with hospital IT networks.

**Backup memory** for easy data collection.

**Data buffer** for paperless data management.

**Maternal and fetal monitoring during transport** with optional battery.

**Transducer plug and play** with automatic screen layout.

**Monitor triplets** with advanced technology (optional).

**Bright and colorful 6.5" (16.51 cm) display** with large numerics and graphics.

**Intuitive touchscreen** for ease of use.

**Printed report after NST** available on fetal heart rate strip automatically or on demand.

**Smart transducers** measure maternal and fetal pulses separately with cross-channel verification for automated coincidence detection.





Battery option facilitates continuous monitoring during maternal transport in healthcare facilities.



New Avalon Toco MP transducers have the ability to monitor the maternal pulse to allow continuous use of the cross-channel verification feature, even when not monitoring maternal SpO<sub>2</sub> or ECG.



Avalon FM series monitors are compatible with the Avalon Cordless Transducer System (CTS) to provide freedom of movement for the mother and continuous monitoring by the care team.

	FM20	FM30
Care stage	Antepartum	Intrapartum
Patient type	Fetal and maternal	Fetal and maternal
Waveforms (not for diagnostic use)	MECG	DECG, MECG
External fetal parameters	US, Toco	US, Toco
Twin capability	Standard	Standard
Triplets capability	Optional	Optional
Internal fetal parameters	N/A	DECG, IUP
Maternal parameters	MECG, NIBP, Pulse	MECG, NIBP, Pulse, SpO <sub>2</sub>
Smart Pulse technology	Standard	Standard
Cross channel verification	Standard	Standard
Fetal movement profile	Standard	Standard
System interface (optional)	Serial, LAN	Serial, LAN
PS/2 interfaces	Optional	Optional
Monitor screen display	6.5" (16.51 cm)	6.5" (16.51 cm)
Touch screen operation	Standard	Standard
NST timer	Standard	Standard
NST trace interpretation (optional)	For up to three FHRs	For up to three FHRs
Data buffer (standard)	Up to one hour	Up to one hour
Battery operating time (optional)	Up to four hours	Up to four hours
Weight	11.2 lbs (5.1 kg)	11.2 lbs (5.1 kg)
Weight with battery option	11.7 lbs (5.3 kg)	11.7 lbs (5.3 kg)

# Avalon FM20 and FM30 fetal monitors fit right in



M2702A Avalon FM20

M2703A Avalon FM30

A wide range of mounting solutions and high-quality Philips medical supplies optimize performance.

Please visit [www.philips.com/AvalonFM20](http://www.philips.com/AvalonFM20)



© 2012 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Printed in The Netherlands  
4522 962 86461 \* JUN 2012